

Instructions for use



Anaftin® Gel

Anaftin® Gel helps in the management of pain caused by small lesions in the mouth, such as those due to aphthous stomatitis, aphthous ulcers, as well as small lesions from braces and ill-fitting dentures. It is also indicated for diffuse aphthous ulcers.

Mode of action:

Anaftin® Gel contains polyvinylpyrrolidone (PVP) and hyaluronic acid, which form a protective film that adheres to the oral cavity, providing a mechanical barrier to the affected areas, thus relieving the pain caused by exposed nerve endings. The hyaluronic acid and Aloe vera in the formulation support the natural healing of damaged tissues.

Instructions for use:

Apply 1 or 2 drops of gel in order to cover the whole mouth ulcer or lesion. Avoid direct contact of the applicator with the lesion. Do not touch the lesion with the tongue for at least 2 minutes so as to allow the protective film to form. Replace the cap after each application. For best results, avoid eating or drinking for at least 1 hour after application.

Use 3 – 4 times a day, or more often if needed.

Ingredients:

Aqua, polyvinylpyrrolidone (PVP), maltodextrin, propylene glycol, PEG-40 hydrogenated castor oil, xanthan gum, potassium sorbate, sodium benzoate, sodium hyaluronate, aroma, benzalkonium chloride, disodium EDTA, sodium saccharin, dipotassium glycyrrhizate, *Aloe barbadensis*

Cautions:

Each Anaftin® Gel tube should be used only by one single patient to avoid cross-contamination of the applicator and gel. The use of Anaftin® Gel is contra-indicated in any patient with a known history of hypersensitivity to any of the ingredients. If you're pregnant ask a doctor before using the product. If you frequently suffer from mouth ulcers, it is recommended that you seek medical advice.

Precautions and warnings:

Do not use after the expiry date, or if the tube appears damaged or is leaking. Keep out of the sight and reach of children. Store at room temperature, away from direct heat. Do not freeze.

Contents:

8 ml

Manufacturer:

Alliance Pharma Srl,
Viale Restelli 5,
20124 Milan,
Italy



Distributed by

BERLIN-CHEMIE AG (MENARINI GROUP),
Glienicke Weg 125,
12489 Berlin,
Germany

Rev 09/2016

engl. / V112939/01
TR1273812 0ALL2FI00044