



ΕΝΔΕΙΞΕΙΣ

ΣΥΣΤΑΤΙΚΑ
AQUA, HYDROGENATED STARCH HYDROLYSATE,
GLYCERIN, PVP, SODIUM CITRATE, POLYSORBATE 80,
PVM/MA COPOLYMER, MALTODEXTRIN, XANTHAN
GUM, PENTYLENE GLYCOL, GLUCONOLACTONE,
SODIUM HYDROXIDE, PEG-40 HYDROGENATED
CASTOR OIL, SODIUM HYALURONATE, SODIUM
BENZOATE, ALOE BARBADENSIS LEAF JUICE,
GLYCERYL CAPRYLATE/CAPRATE, POLYCARBOPHIL,
ALTHAEA OFFICINALIS ROOT EXTRACT, SUCRALOSE,
AROMA, CALCIUM GLUCONATE, POTASSIUM
SORBATE.

ΤΡΟΠΟΣ ΧΡΗΣΗΣ

Για μικρές μεμονωμένες πληγές. Μετά τη σωστή στοματική υγιεινή, εφαρμόστε 1 ή 2 σταγονίους γέλης ώστε να καλυφθεί όλο το έλκος ή η πληγή στο στόμα. Αποφύγετε την άμεση επαφή του απλικατέρ με την πληγή. Εάν είναι δύσκολο να προσεγγίσετε την πληγή, εφαρμόστε τη γέλη με καθαρό δάκτυλο ή γόδα. Μην ακουμπάτε την πληγή με τη γλώσσα για τουλάχιστον 2 λεπτά ώστε να σχηματιστεί το προστατευτικό φιλμ. Κλείνετε το καπάκι μετά από κάθε εφαρμογή. Χρησιμοποιήστε 3-4 φορές την ημέρα ή όσο χρειάζεται.

Χρησιμοποιήστε 3-4 φορές την ημέρα ή όσο χρειάζεται.

ΠΡΟΕΙΔΟΠΟΙΗΣΕΙΣ

ΦΥΛΑΞΗ

Μέγεθος: σωληνάριο 8 ml e

 Biokosmes s.r.l.
Via Besana 10,
20122 Milano ITALY

Αποκλειστικός Διανομέας :
 AMBITAS Healthcare
 Νικηταρά 4 & Ψαρρών, Χαλάνδρι 15232, Ελλάδα
 Τηλ.: +30.210.682.8810
www.ambitashealthcare.com



ALLIANCE
fMBITAS
making healthcare work

PAIN RELIEF OF MOUTH ULCERS AND SMALL ORAL LESION

INDICATION

Aloclair® Plus Bioadhesive Soothing Gel has a mechanical action (barrier effect) indicated for the management of the symptoms caused by aphthous ulcers, oral stomatitis, mouth sores, oral irritations and small oral injuries due to braces and dentures. Its mechanism of action leads to relief of pain by adhering to the mucosal surface of the mouth, soothing oral lesions, reducing pain and discomfort. The gel structure of Aloclair® Plus Bioadhesive Soothing Gel, together with its special applicator, is particularly indicated for clearly visible and identified lesions. The Gel is formulated with ingredients with high film forming power and highly mucoadhesive: it creates a protective film as barrier effect (or a shield) which lasts more than 12 hours. It protects from external agents, minimizing the risk of further irritation. The Gel adheres to mouth ulcers and lesions limiting the contact with external agents such as chemical, mechanical and thermal stimulus. Thanks to its protective action, it helps in the management of the physiological process of healing.

COMPOSITION

AQUA, HYDROGENATED STARCH HYDROLYSATE, GLYCERIN, PVP, SODIUM CITRATE, POLYSORBATE 80, PVM/MA COPOLYMER, MALTODEXTRIN, XANTHAN GUM, PENTYLENE GLYCOL, GLUCONOLACTONE, SODIUM HYDROXIDE, PEG-40 HYDROGENATED CASTOR OIL, SODIUM HYALURONATE, SODIUM BENZOATE, ALOE BARBADENSIS LEAF JUICE, GLYCERYL CAPRYLATE/CAPRATE, POLYCARBOPHIL, ALTHAEA OFFICINALIS ROOT EXTRACT, SUCRALOSE, AROMA, CALCIUM GLUCONATE, POTASSIUM SORBATE.

HOW TO USE

For small single lesions. After a correct oral hygiene, apply 1 or 2 drops of gel so as to cover the whole mouth ulcer or lesion. Avoid direct contact of the applicator with lesion. If the lesion is difficult to reach, use a clean finger or gauze to apply the gel. Do not touch the lesion with the tongue for at least 2 minutes to allow the formation of the protective film. Close the cap after each application. Use 3-4 times a day, or as needed, until the symptoms disappear. After thirty days of use, suspend the

application for at least one week, before restarting the treatment. For best results, do not eat or drink for at least one hour after application.

WARNINGS

Carefully close the cap after each use. There are no known side effects deriving from the use of the product, however it is contra-indicated in patients with a known history of hypersensitivity to any of the ingredients. If the symptoms persist or in case of symptoms worsening, ask for a medical advice. The product is not intended to be ingested, but topically applied in the oral cavity. If the product is swallowed accidentally, no adverse effects are anticipated. Do not use after the expiry date, which is indicated in the pack. Keep out of the sight and reach of children. Do not use if the package is damaged or open. The product should be used in a single patient to avoid cross contamination. Avoid contact with eyes. The product has not been evaluated during pregnancy and breastfeeding. If you are pregnant/breastfeeding, consult your doctor before use. Suitable for children from 4 years of age.

STORAGE

The expiry date refers to the intact product correctly stored. Store in the original pack, in a dry place away from light and heat source. Do not refrigerate the product. Store at room temperature.

Size: 8 ml  tube

Aloclair®, Alliance and medical devices associated are registered trademarks of the Alliance Pharma PLC Group®.

 Biokosmes s.r.l.
Via Besana 10,
20122 Milano ITALY

Exclusive Distributor:
AMBITAS Healthcare
4, Nikitara & Psarron str., 152 32, Chalandri, Hellas
Tel.: +30.210.682.8810
www.ambitashealthcare.com

Symbol		 0477				
Meaning	Legal manufacturer of medical device	CE mark from notified body	Consult the instruction for use for important information	Expiry date	Batch number	Consult the instruction for use for important information